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Table A. Filtration Efficiency Requirements Per U.S., European, and China Standards

Standard	Classification	Filter Efficiency
NIOSH 42 CFR 84 (US)	N95	≥ 95%
NIOSH 42 CFR 84 (US)	N99	≥ 99%
NIOSH 42 CFR 84 (US)	N100	≥ 99.97%
EN 149:2001	FFP1 (filtering facepiece)	≥ 80%
EN 149:2001	FFP2 (filtering facepiece)	≥ 94%
EN 149:2001	FFP3 (filtering facepiece)	≥ 99%
EN 143:2000, EN 140:1999, EN136:1998	P1 (elastomeric facepiece)	≥ 80%
EN 143:2000, EN 140:1999, EN136:1998	P2 (elastomeric facepiece)	≥ 94%
EN 143:2000, EN 140:1999, EN136:1998	P3 (elastomeric facepiece)	≥ 99.95%
GB2626-2006	KN/KP90	≥ 90%
GB2626-2006	KN/KP95	≥ 95%
GB2626-2006	KN/KP100	≥ 99.97%

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	CLOTH FACE COVERINGS	Surgical Masks	N95 Respirators
Testing and Approval	None	Cleared by US Food and Drug Administration (USDA)	Evaluated, tested, and approved by NIOSH as per the requirements in 42 CFR Part 84
Intended Use	Daily use to create a barrier and limit disease spread.	Provides wearer protection against fluids. Protects patient from wearer's respiratory emissions.	Reduces exposure to particles including small particle aerosols and large droplets (only non-oil aerosols).
Face Seal and Fit	Loose	Loose	Tight
Fit Testing Required	No	No	Yes
Respiratory Protection	No	No	Yes (filters 95% of particulates)
Use Limitations	May be re-used with proper care and washing.	Disposable. Must be discarded after each encounter.	Ideally should be discarded after each use. Especially if it becomes wet, dirty, or loose.

NIOSH Particulate Filter Approval Categories

Minimum Filter Efficiency	N series Not resistant to oil	R Series Somewhat resistant to oil	P Series Strongly resistant to oil
95%	N95	R95	P95
99%	N99	R99	P99
100% (99.97%)	N100	R100	P100 (~HEPA)

Just In Time Training

PPE: Respirator Usage and Safety

Understanding the Difference

Surgical Mask

N95 Respirator

Testing and Approval

Cleared by the U.S. Food and Drug Administration (FDA)

Evaluated, tested, and approved by NIOSH as per the requirements in 42 CFR Part 84

Intended Use and Purpose

Fluid resistant and provides the wearer protection against large droplets, splashes, or sprays of bodily or other hazardous fluids. Protects the patient from the wearer's respiratory emissions.

Reduces wearer's exposure to particles including small particle aerosols and large droplets (only non-oil aerosols).

Face Seal Fit

Loose-fitting

Tight-fitting

Fit Testing Requirement

No

Yes

User Seal Check Requirement

No

Yes. Required each time the respirator is donned (put on)

Filtration

Does NOT provide the wearer with a reliable level of protection from inhaling smaller airborne particles and is not considered respiratory protection

Filters out at least 95% of airborne particles including large and small particles

Leakage

Leakage occurs around the edge of the mask when user inhales

When properly fitted and donned, minimal leakage occurs around edges of the respirator when user inhales

Use Limitations

Disposable. Discard after each patient encounter.

Ideally should be discarded after each patient encounter and after aerosol-generating procedures. It should also be discarded when it becomes damaged or deformed; no longer forms an effective seal to the face; becomes wet or visibly dirty; breathing becomes difficult; or if it becomes contaminated with blood, respiratory or nasal secretions, or other bodily fluids from patients.

Centers for Disease Control and Prevention
National Institute for Occupational Safety and Health

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