






UROGEN WELL D-ONE® COLOR TABLE

MYCOPLASMA SECTOR		
 <p>ID ID ID ID <i>M. hominis</i> <i>M.spp</i></p>		
NORMAL IMAGE OF <i>Mycoplasma</i> SECTOR		
MOST COMMON POSITIVE REACTIONS		
UROGEN WELL D-ONE® PLATE PHOTO	INTERPRETATION	SCHEMATIC IMAGE OF REACTIONS
 <p>ID ID ID ID <i>M. hominis</i> <i>M.spp</i></p>	<p>Presumptively <i>Mycoplasma hominis</i></p>	 <p>ID ID ID ID <i>M. hominis</i> <i>M.spp</i></p>
 <p>ID ID ID ID <i>M. hominis</i> <i>M.spp</i></p>	<p>Presumptively <i>Mycoplasma spp.</i></p>	 <p>ID ID ID ID <i>M. hominis</i> <i>M.spp</i></p>

UROGEN WELL D-ONE® COLOR TABLE





UREAPLASMA UREALITICUM/PARVUM SECTOR



U. urealyticum / parvum

NORMAL IMAGE OF *Ureaplasma urealyticum/ parvum* SECTOR

MOST COMMON POSITIVE REACTIONS



UROGEN WELL D-ONE® PLATE PHOTO	INTERPRETATION	SCHEMATIC IMAGE OF REACTIONS
 <p>10⁴ ≥10⁵ ID <i>U. urealyticum / parvum</i></p> <p>GC Growth Control</p>	<p>Presumptively <i>Ureaplasma urealyticum / parvum</i> ≥ 10⁵</p>	 <p><i>U. urealyticum / parvum</i></p> <p>Growth Control</p>
 <p>10⁴ ≥10⁵ ID <i>U. urealyticum / parvum</i></p> <p>GC Growth Control</p>	<p>Presumptively <i>Ureaplasma urealyticum / parvum</i> 10⁴</p>	 <p><i>U. urealyticum / parvum</i></p> <p>Growth Control</p>

Enterococcus spp. SECTOR



NORMAL IMAGE OF *Enterococcus spp.* SECTOR

MOST COMMON POSITIVE REACTIONS

UROGEN WELL D-ONE® PLATE PHOTO	INTERPRETATION	SCHEMATIC IMAGE OF REACTIONS
 <p>ID <i>Enterococci</i></p>	<p>Presumptively <i>Enterococcus spp.</i></p>	 <p><i>Enterococci</i></p>



UROGEN WELL D-ONE® COLOR TABLE

ESCHERICHIA COLI SECTOR



NORMAL IMAGE OF *Escherichia coli* SECTOR

MOST COMMON POSITIVE REACTIONS



UROGEN WELL D-ONE® PLATE PHOTO	INTERPRETATION	SCHEMATIC IMAGE OF REACTIONS
	<p>Presumptively <i>Escherichia coli</i></p>	 <p><i>E.coli</i></p>

Streptococcus agalactiae SECTOR



NORMAL IMAGE OF *Streptococcus agalactiae* SECTOR

MOST COMMON POSITIVE REACTIONS

UROGEN WELL D-ONE® PLATE PHOTO	INTERPRETATION	SCHEMATIC IMAGE OF REACTIONS
	<p>Presumptively <i>Streptococcus agalactiae</i></p>	 <p><i>Strepto B</i></p>

UROGEN WELL D-ONE® COLOR TABLE



STAPHYLOCOCCUS AUREUS SECTOR



ID
S. aureus

NORMAL IMAGE OF STAPHYLOCOCCUS SPP SECTOR.

MOST COMMON POSITIVE REACTIONS



UROGEN WELL D-ONE® PLATE PHOTO	INTERPRETATION	SCHEMATIC IMAGE OF REACTIONS
 <p>ID <i>S. aureus</i></p>	<p>Presumptively <i>S. aureus</i></p>	 <p><i>S. aureus</i></p>

CANDIDA spp. SECTOR



CAN

NORMAL IMAGE OF CANDIDA spp

UROGEN WELL D-ONE® PLATE PHOTO	INTERPRETATION	SCHEMATIC IMAGE OF REACTIONS
 <p>CAN</p>	<p>Presumptively <i>C. albicans</i></p>	 <p>CAN</p>



UROGEN WELL D-ONE® COLOR TABLE

GARDNERELLA VAGINALIS SECTOR



NORMAL IMAGE OF *Gardnerella vaginalis* SECTOR

MOST COMMON POSITIVE REACTIONS

UROGEN WELL D-ONE® PLATE PHOTO	INTERPRETATION	SCHEMATIC IMAGE OF REACTIONS
	Presumptively <i>Gardnerella vaginalis</i>	 <i>Gardnerella vaginalis</i>



NEISSERIA spp. SECTOR






Neisseria spp

NORMAL IMAGE OF *NEISSERIA spp.* SECTOR

MOST COMMON POSITIVE REACTIONS

UROGEN WELL D-ONE® PLATE PHOTO	INTERPRETATION	SCHEMATIC IMAGE OF REACTIONS
 <i>Neisseria spp</i>	Presumptively <i>Neisseria spp.</i>	 <i>Neisseria spp</i>

ANTIBIOTIC SUSCEPTIBILITY TEST SECTOR (*M.hominis* and *U.urealyticum/parvum*)

UROGEN WELL D-ONE® PLATE PHOTO	INTERPRETATION	SCHEMATIC IMAGE OF REACTIONS
	RESISTANCE	
	SUSCEPTIBLE	