

Spokeman Mini-Buddy Manual

Product Overview:

Top Profile



Headset Microphone Port

Side Profile



- 1: Power indicator light: Green when there is sufficient power.
Red when battery is weak.
- 2: On/Amplification/Off knob: Turn clockwise to increase volume.
- 3: Microphone port: Standard 3.5mm mono jack.
- 4: Charging indicator: Red when charging is in progress.
- Green when charge is complete.
- No light when charger is not connected.
- 5: Charging port: Use only the ALDS charger provided.

Specifications:

- * Size Dimensions: 85 (H) x 85 (W) x 35 (D) mm
- * Weight: **175** grams (including battery)
- * Operating Voltage: **6 v**
- * Up to 8 watts of audio power
- * Up to 20dB with neck worn microphone **HM900NI**

Distributed By:

MiniBuddy

Instruction Manual



ALDS
MiniBuddy

The Worlds Smallest
High Powered Voice Amplifier

Spokeman Mini-Buddy Manual

Introduction:

Congratulations on your purchase of an ALDS MiniBuddy

The ALDS MiniBuddy Voice Amplifier has been developed for those with temporary or permanent voice impairment. When using the MiniBuddy a comfortable speech level can be achieved for effective everyday communications while reducing vocal stress and strain.

The ALDS MiniBuddy can boost vocal output up to 6 times without audio feedback when used with the new improved ALDS rear wear headset boom microphone (MB8HM) or the lightweight neck-worn microphone (NM900NI).

- * Featuring environmentally friendly high capacity Nickel Metal Hydride batteries (NiMH) that can be safely recharged overnight.
- * Simple on / amplification / off switch for ease of use.
- * Lanyard, belt mount or shoulder strap convenience.
- * Single power & charging indicator light.
- * MiniBuddy cover protector with beltloop.

The MiniBuddy allows dependable operation from 4 to 8 hours of continuous use depending on amplification requirements. The greater the amplification, the greater the draw on the battery.



MiniBuddy - Warranty Instructions:

Warranty: One Year Limited warranty

ALDS warrants the MiniBuddy personal voice amplifier against defects in materials and workmanship for a period of 12 months from the date of purchase. The NiMH batteries and microphones are warranted for a period of 6 months from date of purchase.

The warranty applies to the original user and the term is established by the original purchase receipt or invoice.

Exclusions:

The warranty does not cover damage to or failure of the unit components due to heat, moisture, shock, misuse, abuse, poor maintenance, improper charging, unusual wear and tear or tampering.

Warranty Repair or Replacement:

In the event of failure, the MiniBuddy voice amplifier, batteries, microphone, or other accessories will be repaired or replaced at the manufacturer's discretion.

Warranty Procedure:

- Contact your dealer for assistance. If you need help locating a dealer in your area, contact ALDS at 1-866-845-2537 for referral or assistance.
- Systems being returned for service should be complete with all components and accessories for review.
- All systems being returned for service or review must be carefully packaged to protect the equipment.
- Include your name, address, day time telephone number, a copy of your purchase receipt or invoice and a written description of the problem.
- Ship by a traceable method to insure arrival and to protect your investment.

Factory warranty / repair depot:



3180-216th. St. Langley, BC V2Z 2E6