



FONA – Every dental solution you need

dental profession, we deliver complete, reliable and accessible is constantly on the move, meeting partners and dentists, so we solutions. Regardless of country or specialisation, every dentist can understand their needs and develop products that support can find a fitting solution in FONA today. As our knowledge and enhance their daily work.

Through decades of experience and deep understanding of the and expertise grow, so does our range. Our international team

Purgenius

SAFE, PURIFIED WATER FOR YOUR TREATMENT





CHAIRSIDE CAD/CAM

TRADITIONAL DENTISTRY





INSTRUMENTS AND LASER





Headquarters DISTRIBUTOR

FONA Dental s.r.o. Ševčenkova 34 85101 Bratislava Slovak Republic

info@fonadental.com www.fonadental.com







Purgenius

Purgenius

MORE THAN PURE WATER

Purgenius is the solution for your safer treatments. It is a unique product, designed for external water filtration, softening and disinfection. The disinfection process is free of added chemicals and does not produce waste water.

COMPLETE PROCESS

- Removing 100 µm particles -
- Removing old chlorine residue and colour -
- Removing calcium -
- Removing 1 µm particles -
- Producing disinfection -
- Killing microorganisms / biofilm

PURE WATER - PROTECTING YOUR HEALTH AND INVESTMENT

Removing microorganisms and biofilm protects the health of you and your patients. Eliminating particles and calcium prevents sedimentation and thus protects your instruments and treatment unit.



Purgenius Solution for one or two treatment units.



Purgenius with Booster Water solution for up to six treatment units.

FEATURES AND BENEFITS

Universal Disinfection System

Purgenius can be connected to any FONA treatment unit. As an external device, it is completely independent from the treatment center water unit size and suction system used. With flexible installation, it can be connected chair side or outside of the treatment room.

Chemical free – Harmless Solution

The unique use of free chlorine, produced by Purgenius from tap water as disinfectant, is both ecological and safe for the patient and user.

Sound Investment for Your Health

Not only patients but also dentists are exposed to microrganisms found in instrument cooling water. Reducing microorganisms protects the health of both, while at the same time removing calcium and other unwanted minerals, protecting your instruments and treatment unit from issues caused by sediments.

TECHNICAL SPECIFICATION

DIMENSIONS	280
ELECTRICAL POWER	110
INCOMING WATER PRESSURE	2-6
WATER FILTER CAPACITY	ma
WATER FLOW OUT	1,5
OUTLET WATER PRESSURE	3 b

SCOPE OF SUPPLY

MAIN UNIT	Pur
POWER CABLE	EUt
DOCUMENTATION	Inst
ACCESSORY	Adju
STARTER KIT	Wat



80 mm (l), 166 mm (w), 380 mm (h)

10-230 V, 50/60 Hz, 200 W max (<10 W standby)

-6 bar

ax 2.200 l

i-2 l/min

oar (adjustable)

rgenius

type

structions for use, Installation and service form, test form

justment pin and sealing cup

ter Filter , Bio reduction filter, chlorine test strips