

PELICAN WHEELCHAIR

BOOT SLIDER



- People often have difficulty getting a wheelchair into the boot of a car, and usually resort to having to lift the wheelchair. The Wheelchair Boot Slider can help make this job easier.
- The Boot Slider consists of three parts.
- The Boot Blanket is the padded section that protects the rear of the car. The padded section goes inside the boot to give a slippery surface for the wheelchair to slide on, when putting it into or taking it out of the car.
- The Slider Platform is what the wheelchair is rolled into. It has a cover to hold the wheelchair in position.
- Place the Ground Protector on the ground under the Slider Platform to protect the Slider Platform from getting dirty, wet, or deteriorating on rough road surfaces. It also keeps your clothes clean when leaning the Wheelchair Boot Slider up against your thigh. The Ground Protector has a webbing strap attaching it to the Boot Blanket. This is a reminder not to leave it behind before driving away!



Warranty
1 Year

Compliance
No applicable compliance standard

Code	Description	Weight
580K	Pelican Wheelchair Boot Slider	6kg