



Northern Lite is a lightweight, compact, folding wheelchair that is easy to adjust for optimal performance and functionality. With a center of gravity that can be adjusted in 1/4" increments, plus a vertically adjustable axle plate, this is one of the most customizable folding wheelchairs on the market today.



### **Features**

Adjustable camber angle on rear wheels Center of gravity adjustments in ¼" increments Composite slip-up footplates Caster options 4", 6" or 8" Adjustable tension back upholstery Seat depths 16" - 18" Removable, height adjustable T-style armrest Ultra-light 70° pin style footplates Front seat to floor height 14" - 19" Multi position caster fork Seat widths 16" - 20" Back heights 12" - 20" Seat upholstery bottom storage pocket



perfect for today's 'on the go' user.

#### **REAR AXLE PLATE**

The Axle plate allows the user to adjust their center of gravity, it also allows many different seats to floor heights.



#### **QUICK FLIP ANTI TIPPERS**

Easy to use quick flip and out of the way when required, while providing safety that is fall proof prevention system designed to protect you and your loved ones from harm.



#### **DOUBLE CROSS BRACE**

The Northern Lite uses Double Cross brace construction. This feature improves strength and rigidity of Northern Lite.



#### WHEEL LOCK

Our lock system is designed to keep your wheelchair in 'park' mode when performing daily tasks that involve upper body movement.









**Seat Width** 10" - 20", 22" - 24"



Caster Size 4", 5", 6", 7" Urethane



Rear Wheel 20", 22" or 24" Airless or Pneumatic



**Seat Depth** 12" - 20"



Back Style Slip on Back Upolstery



Arm Style Adj. Height T-Style Gel or Vinyl Adj. Height Flip Back Gel or Vinyl



Front Rigging
Dual Swingaway, Dual Elevating
Legrest, or Elevating Legrest



Seat to Floor 13" - 16.5" Small Frame 17" - 19.5" Large Frame



Frame Material Lightweight Aluminum



Product Weight 32 lbs. (14.5 kg)



Overall Width + 9.5"



Weight Capacity 250 lbs. (113kg)

ADP Code WM0001188 HCPCS Code K0005



# Warranty:

Future Mobility Healthcare Inc. warrants the Northern Lite wheelchair to be free from defects in materials and workmanship for two (2) years from the date of purchase. The frame has a lifetime warranty and the upholstery for 90 days, upon normal usage by the original purchaser. If within this warranty period the product shall be proven to be defective, such product shall be repaired or replaced, at Future Mobility Healthcare Inc. discretion.





6750 Professional Court, Mississauga, ON L4V1X6

1-888-737-4011

www.futuremobility.ca

## USA

1 Buffalo River Place, Buffalo, NY 14210

1-855-938-8873

www.futuremobility.com