INVACARE SOFTFORM PREMIER **ACTIVE 2RX HYBRID MATTRESS**





- High/Very High Risk. The Softform Premier Active 2 retains properties as a high risk static mattress but should clinical judgement determine that the patient requires stepping up to a dynamic surface.
- A discrete pump can be fitted to the air insert transforming the mattress into a dynamic surface and delivering additional levels of pressure relief.
- Using the pump allows stepping up and down between reactive and active alternating systems with minimum patient disruption
- Inflation of air cells alternates on a 10 minute cycle promoting tissue recovery and pressure distribution
- In the event of a power failure, the patient is still resting on a high risk pressure reducing surface



Warrantv 2 Years

Compliance

EN 597-1, EN 597-2 and BS 7177 Crib 5 | EN 60601-1

Code	Description	Width	Length	Height	S.W.L
PKG-ACTIVE2RX	Invacare Softform Premier Active 2 RX Mattress with Pump - Hybrid	880mm	1970mm	152mm	247.6kg

Phone: 07 4208 8166 Email: sales@peakmobilityandrehab.com.au