



ORION II TILT-IN-SPACE

Future Mobility's Orion II is one of the most innovative, popular tilt wheelchair on the market today. It was developed for active individuals who need to redistribute their body weight to achieve optimal position.



A STANDARD WEIGHT LIMIT OF 350 LBS.




STANDARD WEIGHT
LIMIT 350 LBS

NO ADDITIONAL COST

Our Orion II Wheelchairs are up to 350 LBS standard, with
NO Additional Cost.





13.25"
**SEAT-TO-FLOOR
HEIGHT**

ANTI-TIPPER

Safety is essential. That's why our Anti-Tipper fall prevention system is sold standard with all Future Mobility wheelchairs. The system is designed to keep your loved ones safe from harm.



ARMREST

Future Mobility's wheelchairs are sold with a variety of armrest options. This ensures the end user has the most comfort and support possible when using our wheelchair. Our Armrests are available in Flip-up, height adjustable, and/or fixed style.



POSI-LOCK

Located underneath the chair, the Posi-Lock allows for a gentle and slow tilt of the end user. It is lockable at different angles.

SWINGAWAY FRONT RIGGINGS

Swingaway Front Riggings give the end user the option of swing-out or dual swing-in and swing-out. Ideal for transporting and seating the end user.




SPECIFICATIONS:

ORION II TILT-IN SPACE



SEAT TO FLOOR

13 1/4" - 20"




REAR WHEEL

12", 20", 22" OR 24"
COMPOSITE URETHANE



SEAT WIDTH

15" - 24"



SEAT DEPTH

15" - 24"



OVERALL WIDTH

SEAT WIDTH +9.5"



CASTER SIZE

4", 5", 6", 7" AND 8" URETHANE



FRONT RIGGING

PINSTYLE, SWINGAWAY
DUAL SWINGAWAY
OR DUAL ELEVATING LEGRESTS



PRODUCT WEIGHT

65 LBS (29KG)



ADJUSTABLE BACK ANGLE

-7°, 0°, 7°, 14°, 21°



BACK STYLE

ORON BACK



ARM STYLE

ADJ. HEIGHT T-STYLE
ADJ. HEIGHT PIN RELEASE T-STYLE
ADJ. HEIGHT FLIP BACK



WEIGHT CAPACITY

350 LBS (159 KG)



ADP
Ontario

CODE

WMA5FM005

WARRANTY:

The Capella, Capella45, Orion II, Orion III and the Super Chair for Kids is warranted to be free from defects in materials and workmanship for a period of one (1) year and the Orion II 500 for period of one (1) year from the date of purchase. The frame has a lifetime warranty and the upholstery for 90 days, upon normal usage by the original purchaser. Power tilt batteries are warranted for six (6) months while the control box, battery charger, actuator and hand control are warranted for one (1) year. If within this warranty period the product shall be proven to be defective, such product shall be repaired or replaced, at Future Mobility Healthcare Inc. discretion.



CANADA

6750 Professional Court, Mississauga,
Ontario, L4V 1X6



1-888-737-4011



www.futuremobility.ca



USA



1 Buffalo River Place, Buffalo,
NY, 14210



1-855-938-8873



www.futuremobility.com