ASPIRE OVER TOILET AID







aspire

- Raised armrests Assist the user transferring on and off the Over Toilet Aid
- Additional height increments on one leg Accommodates uneven surfaces
- Height adjustable legs Easy push button height adjustment for changing needs
- Comfortable, contoured seat with clip on attachment For easy assembly and cleaning. Padded seat available, sold separately (BTT145110)
- All frames are individually treated for corrosion resistance
- Removable, moulded splash guard provided Easily removed for transport and cleaning
- Durable non-slip feet Provide excellent grip and stability in wet conditions. Stainless steel washer within enhances durability
- Available in Aluminium or Treated Steel

Warranty Compliance

5 Years Corrosion 1 Year Components AS/NZS3973.7:2009

Code	Description	Overall Height	Overall Width	Overall Depth	Seat Depth	Seat Width	Between Arms	Seat to Floor Height	Product Weight	S.W.L
BTT145000	Over Toilet Aid - Steel	700 - 880 mm	512 mm	625 mm	435 mm	370 mm	460 mm	450 - 630 mm	4.2 kg	175 kg
BTT145005	Over Toilet Aid - Aluminium	700 - 880 mm	512 mm	625 mm	435 mm	370 mm	460 mm	450 - 630 mm	3.4 kg	175 kg
BTT145110	10 Over Toilet Aid - Padded Seat - Aspire - Grey									
BTT145105 Over Toilet Aid - Plastic Splash Guard only										

Phone: 07 4208 8166 Address: Shop 8C, 58-62 Norman Street, Gordonvale Qld 4865

Email: sales@peakmobilityandrehab.com.au