



# Clinical Engineering Services

## Lift-Assisted Stirrups



Side View

### BD100

#### Lift-Assisted Stirrups

Travels -30° to +120°, from +90° to abduction & Adduction Done at Handle  
Includes Boot Liners & Side Rail Clamps  
500 lb. Weight Capacity  
Pair



### DSC-800-0342-PL

#### Great White Platinum Stirrups

Range of +80° Motion: -30° to abduction -9° to +25° & Adduction  
Includes Premium Boot Liners  
600 lb. Weight Capacity  
Optional: Clamps Sold Separately  
Pair



### DSC-800-0342-PR (No Fin)

#### Great White Premium Stirrups

Range of +80° Motion: -30° to abduction -9° to +25° & Adduction  
Includes Premium Boot Liners  
400 lb. Weight Capacity  
Optional: Clamps Sold Separately  
Pair



### DSC-800-0342-R

#### Great White Robotic Stirrups

Range of Motion: 65° low to 45°  
High Abduction -9° to +25° & Adduction  
Includes Premium Boot Liners  
600 lb. Weight Capacity  
Optional: Clamps Sold Separately  
Pair



### KTEK3100 (With Fin) / KTEK5100 (With Fin)

#### Lift-Assisted Stirrups

No Velcro® brand fasteners  
Range of +85° Motion: -55° to abduction +25° & Adduction -9°  
KTEK3100: 350 lb. / KTEK5100: 500 lb.  
Includes Boot Liners  
Optional: Clamps Sold Separately  
Pair

[OrderCES@CESsurgical.com](mailto:OrderCES@CESsurgical.com)

Clinical Engineering Services

540-539-3052

[CESsurgical.com](http://CESsurgical.com)