ASPIRE OVER TOILET AID - WIDE









Email: sales@peakmobilityandrehab.com.au

- Raised armrests Assist the user transferring on and off the Over Toilet Aid
- Additional height increments on one leg Accommodates uneven surfaces
- Height adjustable legs Easy push button height adjustment for changing needs
- Comfortable, contoured seat with clip on attachment For easy assembly and cleaning. Padded seat available, sold separately (BTT145110)
- All frames are individually treated for corrosion resistance
- Removable, moulded splash guard provided Easily removed for transport and cleaning.
- Durable non-slip feet Provide excellent grip and stability in wet conditions. Stainless steel washer within enhances durability
- Available in Aluminium or Treated Steel SWL175kgs

Warranty Compliance

5 Years Corrosion 1 Year Components AS/NZS3973.7:2009

Code	Description	Overall	Overall	Overall	Seat	Seat	Between	Seat to Floor	Product
		Height	Width	Depth	Depth	Width	Arms	Height	Weight
BTT145010	Over Toilet Aid - Wide - Steel	700 - 880 mm	570 mm	625 mm	435 mm	370 mm	525 mm	450 - 630 mm	4.2 kg
BTT145015	Over Toilet Aid - Wide - Aluminium	700 - 880 mm	570 mm	625 mm	435 mm	370 mm	525 mm	450 - 630 mm	3.4 kg
BTT145110	Over Toilet Aid - Padded Seat - Asp								
BTT145105	Over Toilet Aid - Plastic Splash Gua								

Phone: 07 4208 8166 **Address:** Shop 8C, 58-62 N

Website: www.peakmobilityandrehab.com.au